

MA-UDS01

Full Digital
Black and White Ultrasound
Diagnostic System





Unmatched Performance
and Precision


Versatile Applications for Comprehensive Care

This device model extends its capabilities across various medical fields, ensuring accurate diagnostics and improving patient care.

Routine Examinations: abdomen and small parts



Specialized Fields: urology, obstetrics, pediatrics



CONVENIENT AND COMPACT DESIGN




15-inch

high-definition LED screen:

180° opening and closing angle range adjustment

Lightweight Body and Simple Design

This device model boasts a lightweight body and a minimalist design for enhanced mobility and user-friendly operation.



Advanced Interfaces and Controls



· 2 probe ports



· A classic operation keyboard



· 1 video interface/ 1 VGA interface/
2 USB interfaces

IMPROVED PERFORMANCE LEADS TO HIGH-QUALITY IMAGING




All-round upgraded chip,
better performance


The new ARM architecture,
improvement of the speed
and capabilities of processing



RICH FUNCTIONS BRING SIMPLIFIER OPERATION




Pseudo color: 3 colors




Pseudo color: 7 colors




Multiple system languages:
CH, EN, ES, FR, RU, PT



Gravel positioning function:
displaying the midline, one-click depth measurement




Puncture guidance function




Rich obstetrics measurement package:
supporting pre-setting function.


CLINICAL IMAGING



B-Mode




B/B-Mode









4B-Mode



M-Mode



B/M-Mode



OPTIONAL PROBES



Convex probe

Abdomen, obstetrics,
gynecology



Linear probe

Vascular,
musculoskeletal



Trans-vaginal probe

Obstetrics, gynecology,
urology

